

Poster Communication 51

GC-MS Terpene Profiling in Cannabis extracts: Method Optimization and Analytical Validation

Albertina M. Capita^{1,2}, **Ricardo J. Dinis-Oliveira**^{1,2,3,4}, **Carlos J. A. Ribeiro**^{1,2,5} and **Eduarda M. P. Silva**^{1,2,*}

¹ Associate Laboratory i4HB - Institute for Health and Bioeconomy, University Institute of Health Sciences - CESPU, 4585-116 Gandra, Portugal.

² UCIBIO - Applied Molecular Biosciences Unit, Translational Toxicology Research Laboratory, University Institute of Health Sciences (1H-TOXRUN, IUCS-CESPU), 4585-116 Gandra, Portugal.

³ Department of Public Health and Forensic Sciences, and Medical Education, Faculty of Medicine, University of Porto, 4200-319 Porto, Portugal.

⁴ FOREN—Forensic Science Experts, 1400-136 Lisbon, Portugal.

⁵ Avextra Portugal SA, 1070-060 Lisboa, Portugal.

* Correspondence: eduarda.silva@iucs.cespu.pt

Abstract

Background: The importance of Cannabis sp. in modern medicine stems from its complex phytochemical profile, where terpenes may play a role through the entourage effect. This entourage effect proposes that Δ^9 -tetrahydrocannabinol (Δ^9 -THC) interacts with other compounds in the cannabis plant in ways that meaningfully modify its effects. For instance, Spindle et al. reported that D-limonene reduced the acute anxiogenic effects of Δ^9 -THC [1]. To ensure the standardization of cannabis-based products, GC-MS is essential by enabling accurate quantification of these compounds in complex plant matrices [2]. **Objective:** This study aims to optimize and validate a GC-MS method for the characterization of major terpenes in cannabis leaves and flowers. **Methods:** Two cultivars of cannabis flowers were used: one rich in CBD (Blue Cheese) and the other rich in Δ^9 -THC (Z-Face). In addition, leaves from a Δ^9 -THC-rich cannabis cultivar provided by Avextra were also analyzed. Samples were pulverized and extracted with ethyl acetate for 10 minutes using a Retsch 400 ball mill. GC-MS analysis was performed using an SH-Rxi-5ms column under a temperature gradient. Terpenes were identified in SIM mode based on characteristic ions [3]. The method was validated according to ICH guidelines. **Results:** The method showed high precision and accuracy, demonstrating linearity ($R^2 > 0.999$) for all ten analyzed terpenes. The three most abundant terpenes identified in the leaves were β -caryophyllene (0.014%), α -bisabolol (0.011%), and α -humulene (0.09%), expressed as percentage (w/w). In the Blue Cheese cultivar, the dominant terpenes were α -bisabolol (0.174%), β -caryophyllene (0.108%), and α -pinene (0.086%). The predominant terpenes in the Z-Face cultivar included β -caryophyllene (0.445%), α -humulene (0.252%), and limonene (0.227%). Terpinolene was the only terpene not detected in all three samples. **Conclusion:** A validated GC-MS method was established for terpene characterization in cannabis extracts. Higher terpene levels were found in flowers compared to leaves, with distinct profiles observed between the CBD- and THC-rich varieties.

Keywords: Cannabis; GC-MS; Terpenes

Acknowledgments/Funding

This work was funded by Avextra Portugal SA. The authors also acknowledge financial support through the annual funding of 1H-TOXRUN of the University Institute of Health Sciences (IUCS-CESPU).

References

1. Spindle, T.R. et al. Vaporized D-Limonene Selectively Mitigates the Acute Anxiogenic Effects of Δ^9 -Tetrahydrocannabinol in Healthy Adults Who Intermittently Use Cannabis. *Drug Alcohol Depend.* **2024**, *257*, 111267, doi: 10.1016/j.drugalcdep.2024.111267.
2. Francisco, V.P. et al. Development of GC-MS Method for the Quantification of Cannabis Terpenes: Application to Five Commercial Varieties. *J Chromatogr B* **2024**, *1247*, 124316, doi.org/10.1016/j.jchromb.2024.124316.
3. Capita, A.M. et al. Analysis of Terpenes in Cannabis by GC-MS: Method Development and its Application to *Cannabis sativa* L. extracts, *Sci Lett*, **2025**, *1* (Sup 1), doi: 10.48797/sl.2025.390.



In *Scientific Letters*, articles are published under a CC-BY license (Creative Commons Attribution 4.0 International License at <https://creativecommons.org/licenses/by/4.0/>), the most open license available. The users can share (copy and redistribute the material in any medium or format) and adapt (remix, transform, and build upon the material for any purpose, even commercially), as long as they give appropriate credit, provide a link to the license, and indicate if changes were made (read the full text of the license terms and conditions of use at <https://creativecommons.org/licenses/by/4.0/legalcode>).